

Fig2-12-3 Method of annotation selection on mRNA variation viewer

Selected annotation is shown

chr1+5241 [Help](#)

Onepass annotation setting window


Color setting Library Seq. ID ☐ Shown ☒ Not shown Abscissa scale 1 [Apply](#)

cDNA annotation setting window

Annotation ☒ Pfam ☒ PROSITE ☒ PSORT ☒ SignalP ☒ SOSUI Abscissa scale 1 [Apply](#)

BC093683.1

Selected all annotations



Selected annotation is shown

chr1+5241 [Help](#)

Onepass annotation setting window


Color setting Library Seq. ID ☐ Shown ☒ Not shown Abscissa scale 1 [Apply](#)

cDNA annotation setting window

Annotation ☐ Pfam ☐ PROSITE ☒ PSORT ☒ SignalP ☒ SOSUI Abscissa scale 1 [Apply](#)

BC093683.1

Selected annotations



Annotation is not shown.

chr1+5241 [Help](#)

Onepass annotation setting window

Color setting Library Seq. ID ☐ Shown ☒ Not shown Abscissa scale 1 [Apply](#)

cDNA annotation setting window

Annotation ☐ Pfam ☐ PROSITE ☐ PSORT ☐ SignalP ☐ SOSUI Abscissa scale 1 [Apply](#)

BC093683.1

Not selected

